




The catalog of products for the visually impaired and hard of hearing
Customer Service: 800-468-4789 9am-5pm Eastern Monday-Friday

Your Minder Talking Pill Reminder Clock LS&S #101068



Buttons on top of clock from left to right (display facing you):

Talk (green) – This button is used to speak the current time.

Alarm check (blue) – This button is used to check what alarms are set.

Light (yellow) – This button turns light on for the display

Alarm Acknowledged (red) - When the alarm is sounding, stops the alarm and tells when next alarm is due.

Buttons on back of clock:

Top slide switch - Alarm selection – This switch is used to select each of the individual alarms during all modes

Left hand side, top to bottom:



The catalog of products for the visually impaired and hard of hearing
Customer Service: 800-468-4789 9am-5pm Eastern Monday-Friday

Alarm off/beep/voice (slide switch) – This switch is used to turn the alarm off, or turn it on to beep or play the pre-recorded voice message. All the way to the left is off, the middle is Beep and all the way to the right is the voice message.

Mode Selection (slide switch): This is used to set the time and the alarm and to “lock” the clock in normal time run mode. All the way to the left is Time Set, the middle is Lock/Run Mode and all the way to the right is Alarm Set.

Set/Advance (push button): You press this button to advance through each step of the setting procedures

"-" and "+" (push buttons): These buttons are used to move selection backward or forward

Right hand side buttons, top to bottom:

Record (red push button): When pressed and held, used to record selected alarm message

Erase (push button): When pressed and held, erases recording for selected alarm

Play (push button): Plays current recorded message for the selected alarm

There is also the ac adaptor jack underneath these buttons.

Setting the current time/date:

Slide Mode Selection switch to time set (all the way to the left). Press

"set/advance" button once to enter Hour Set Mode. You will hear “Hour Set” and current hour setting. Use "-" and "+" buttons to select the desired hour. The clock will announce the hour as you press the buttons.

Once the hour has been set, press the “Set/Advance” button once. You will hear “Minute set” and the current minute setting. Use “-“ and “+” buttons to select



The catalog of products for the visually impaired and hard of hearing
Customer Service: 800-468-4789 9am-5pm Eastern Monday-Friday

the desired minutes. The clock will announce the minute as you press the buttons.

Once the minutes have been set, press the “Set/Advance” button once. You will hear “Month Set” and the current month setting. Use the “-” and “+” buttons to select the desired month. The clock will announce the months as you press the buttons.

Once the month has been set, press the “Set/Advance” button once. You will hear “Day Set” and the current day setting. Use the “-” and “+” buttons to select the desired day. The clock will announce the day numbers as you press the buttons.

Once the day has been set, press the green talk button on the top of the clock to confirm the correct time and date has been set.

Once you have confirmed the correct time and date, slide the Mode Selection switch back to “lock/run” (in the middle).

Setting alarm times:

Slide "Alarm Selection" switch to the desired alarm that you would like to set. All the way to the left is alarm 1. There are 6 alarms that you could set moving the slide switch one click to the right each time.

Slide "Mode Selection" switch to alarm set (all the way to the right). Press the “Set/Advance” button once. You will hear “The (1st, 2nd, 3rd etc) Alarm, hour set” and the current hour setting. Use the “+” and “-” buttons to select the desired hour. The clock will announce the hours as you press the buttons.



The catalog of products for the visually impaired and hard of hearing
Customer Service: 800-468-4789 9am-5pm Eastern Monday-Friday

Once the desired hour has been set, press the “Set/Advance” button once. You will hear “The (1st, 2nd, 3rd etc) Alarm, minute set” and the current minute setting. Use the “+” and “-” buttons to select the desired minute. The clock will announce the minutes as you press the buttons.

Once the desired minute has been set, press the “Set/Advance” button once. You will hear “The (1st, 2nd, 3rd etc) Alarm, (off or on)”. Press either the “+” or “-” button to toggle between on and off.

Once the alarm has been set, you may either move the Alarm Select slide switch to the next alarm and repeat the above steps to set the rest of the alarms, or you may turn the “Mode Selection” slide switch to the Lock/Run position (in the middle). The clock will now be in normal run time mode.

When the alarm is going off, press the red “Alarm Acknowledged” red button (farthest to the right if the display is facing you) to stop the alarm from sounding.

Setting alarm sound:

This clock can play a pre-recorded message as the alarm sound or the alarm can go off as a “beep”.

To set this sound, you will need to slide the Alarm slide switch on the back of the clock to “beep” (in the middle) or “voice” (all the way to the right).



The catalog of products for the visually impaired and hard of hearing
Customer Service: 800-468-4789 9am-5pm Eastern Monday-Friday

Recording your own alarm:

Using the Alarm Select slide switch, select the alarm you would like to record (1-6). You can record your own alarm message by pressing and holding the red “Record” button on the right hand side of the back of the clock. You will hear a “beep” sound indicating that the clock is recording. Speak into the “Mic” speaker that is located directly above the record button. You could say something like: “It is 10AM, remember to take your asprin” or whatever message you would like to play as the alarm goes off. Once you are done recording, release the red “Record” button. You will hear another beep, indicating that the clock has stopped recording.

To hear what you have recorded, press the “Play” button, which is below the “Record” button to the right. You may erase what you have recorded if you are not satisfied or want to change the message by pressing the “Erase” button, which is located below the “Record” button to the left.

You can record a different message for each alarm.

You will want to make sure that the Alarm Sound slide switch is all the way to the right for Voice.

Alarm check:

At any time, you may check to see what time the alarms are set for and if they are on or off. To do this, you would press the blue “Alarm Check” button. This is the second button from the left if the display is facing you. Press the button once to hear the 1st alarm time, twice to hear the 2nd alarm time and so on. This is a good check to do after you have set the alarms or changed the times.



The catalog of products for the visually impaired and hard of hearing
Customer Service: 800-468-4789 9am-5pm Eastern Monday-Friday

You may also press the red “Alarm Acknowledged” button (farthest to the right if the display is facing you) to hear when the next alarm will be going off. It will announce “The next reminder will be the (1st, 2nd, 3rd etc) alarm” and then the time that the alarm is set for.

If you press the “Alarm Acknowledged” button and nothing is announced, none of the alarms have been turned on.

Recording your own alarm messages:

Slide the "Alarm Selection" switch to the alarm you wish to set. (1-6)

Press and hold "record" button when you are ready to record. A 'beep' will notify you that the recording has begun.

- Speak clearly into the "mic" on the back of the clock. (directly above the record button)
- Release "Record" button once you have completed your message. A 'beep' will notify the end of the recording.
- Press "Play" to listen to your recorded message.
- You can re-record or record over the message if it is not to your liking, or you would like to change.
- To stop the alarm while it is going off, press the red Alarm Acknowledged button (far right button on top of clock)