

Settings:

Making Adjustments Press the **MENU** button to display the on screen menu:

Please Note:Changes can only be made to the line when highlighted in brackets and Red font. Use the **UP** or **DOWN** button to highlight the line to apply changes.

Use the **LEFT** or **RIGHT** buttons to Red font the text for adjustment.

Set Time	< 04:47 PM >
Set Date	08-12-2017
Time Mode	12Hours
Date Mode	Day-Month-Year
Language	English
Brightness	on
Alarm Status	
Volume	4
Set Chime	...

Line 1- Set Time:

Press **"OK"** to enter into the time setting. While the digit is Red font the hour<**04**:47>,Press UP or DOWN to change the time.when you get the correct digit,Press **"OK"** to save.Then press **"OK"** and LEFT or RIGHT to highlight the minutes <04:**47**>.Press UP or DOWN to change the time.when you get the correct digit,Press **"OK"** to save(Press **MENU** to exit)

Line2 - Set Date:

Select "Set Date "so that it is highlight in red.Press "OK "to enter into the date setting Font the Day:<**08**-01-2018>or<**01**-08-2018>,use UP or DOWN to adjust date. when you get the correct digit,Press **"OK"** to save, Then press **"OK"** and LEFT or RIGHT to highlight Month/Year <08-**01**-2018>Press UP or DOWN to change the date.when you get the correct digit,Press **"OK"** to save(Press **MENU** to exit)

Line3 - Time Mode:12/24HourTime

Press "OK" to Red font <12 Hours>.LEFT or RIGHT to Select Time Mode between 12-Hours and 24-Hours

After the modifications are completed, Press "OK" to save(Press **MENU** to exit)

Line4 - Date Mode:<Day-Month-Year>or<Month-Day-Year>

Press "OK" to Red font.LEFT or RIGHT to Select Date Mode .

After the modifications are completed, Press "OK" to save(Press **MENU** to exit)

Line5 - Language

Press "OK" to Red font.Then Press LEFT or RIGHT to Select 8 language.

After the modifications are completed, Press "OK" to save(Press **MENU** to exit)

Line6 - Automatic Brightness:

The default is turn on the Auto Brightness, Press "OK" to Red font then Turn to manual adjustment brightness from Level 1 to Level 5 and turn **on** the Auto Brightness .LEFT or RIGHT Select brightness

When you get the brightness you want, Press "OK" to save(Press **MENU** to exit)

Line7 - Alarm Status:

On Line 7 , Press "OK" to enter the alarm settings -- Press "OK" to Add Timer-- Press "OK" to enter Time 1 (Timer1 to Timer 12) alarm setting -- Press "OK" to set Alarm Time, Date, Frequency and Tag

(Tag:Meal, Medicine, Clinic, Shower, Gardening, Cleaning, Appointment and Walk dog)

Note:You need to turn on the **start timer <On>** then the alarm will sound

After the modifications are completed, Press "OK" to save(Press **MENU** to exit)

Line8 - Volume

Press "OK" enter Volume adjustment, from 1--5 The sound will be bigger and bigger

When you select the volume, Press "OK" to save(Press **MENU** to exit)

Line9 - Set Chime

Press **"OK"** enter Hourly time setting, Select the time you want to on-time alarm , press the "OK" button to confirm the selection, and when the number turns red, the selection is successful. (Press **MENU** to exit)

Digital Photo Frame Features:

When you are in the clock interface, after inserting the SD card, press **"OK"** to enter the playback image and video.

support video format: mp4, avi, rmvb.support picture format: Jpg, jpeg

Tip: Please use genuine SD Card and USB Flash Drive,

Maximum support capacity: 64GB