

# BigEye 310 / 518 / 3518

## Assembly Instructions

BigEye Lamp, A Division of Mattingly Low Vision, El Paso, Texas • 888.642.0842



Two lightbulbs and one 5 inch diameter lens are packed in all cartons. LED bulbs are also available for purchase.

### Base Assembly:

1. Stand Base on edge.
2. Insert bolt through hole in the bottom of base.
3. Slide lock-washer over bolt.
4. Align pole with the bolt and thread together.

### Model 310 Height Adjustment:

1. Note: Auto-clutch assembly built in for height adjustment.
2. Loosen the clutch assembly by turning to the left.
3. Set desired height.
4. Turn clutch assembly to the right to tighten.

### Model 518 / 3518 Height Adjustment:

1. Simply pull up or push down to desired height.

### Lens Assembly:

- Note: Take care not to scratch the lens.
- Remove the lens from its carton.
- Remove the bolt and wing-nut from the lens support tongue / connector.
- Line up the holes on the lens and support tongue.
- Insert the bolt from the top and screw on the wingnut (Do not over-tighten).

### Lens Care:

- Only use a lens cloth to clean the lens.
- Do not use window cleaner. Instead, use warm soapy water if necessary.
- The lens is replaceable but is not covered under warranty.

### Proper Use:

- Position your lamp so that you are looking through the lens.
- Adjust the lamp neck up or down as needed to bring the subject matter into sharp focus, then turn the lamp on.
- Take care that the shade and light bulb are not touching fabric, paper, etc.

### PRECAUTIONS:

1. Fire Risk: The lens is a powerful magnifier that can magnify the sun's rays. Do not leave in direct sunlight.
2. Use the lamp's shade for positioning. Do not adjust by twisting the lens. This could break the lens.
3. The lamp will be warm to the touch while in use.
  - a) Give it ventilation space while in use.
  - b) Don't cover the lampshade while in use or warm.
  - c) With prolonged use, the heat from the bulb may produce an odor. This is normal.
  - d) If uncomfortable with the warmth, try the LED bulb. (sold separately)

# BigEye 310 / 518 / 3518

## Warranty and Service

BigEye Lamp, A Division of Mattingly Low Vision, El Paso, Texas 888.642.0842

**All Big Eye products are warranted against manufacturing defects for a period of one (1) year from the date of purchase.**

- Lenses and bulbs are not under warranty.
- Warranty will be void if malfunction results from damage or unreasonable use.

### To Return Lamp for Service:

- Call for repair information:

**BigEye Lamp, Inc**  
A Division of Mattingly Low Vision, Inc  
135 Rio West Drive, Bldg B  
El Paso, TX 79932  
www.mattinglylowvision.com  
888.642.0842

- Prior to any returns, you must contact us. 888-642-0842
- Remove and do not return the lens, lens hardware, and lightbulb from the lamp, or the base and hardware from floor models.
- Enclose a short note with the following information:

- ⇒ Your Name.
- ⇒ Your Complete Address.
- ⇒ Your Telephone Number.
- ⇒ A brief description of the problem.

### Big Eye Products:

Model:	Description:
010	Table Lamp with magnifier
310	Floor Lamp with magnifier
L-465	Power Booster Lens
L-475	Super Booster Lens
L-575	Replacement Lens
PL-40	Replacement Bulb
M14	Free standing magnifier
1418	Deluxe two (2) arm Table Lamp
518	Deluxe Floor Lamp
3518	Deluxe two (2) arm Floor Lamp
LBigEye Rolling Base	4 Wheel Caster Rolling Base
LampBulbLED4k	4000 Kelvin LED Bulb for BigEye Lamp ( 1 year warranty )
LampBulbLED3k	3000 Kelvin LED Bulb for BigEye Lamp ( 1 year warranty )